empelvic®

Instructions for use • Mode d'emploi •
Gebrauchsanweisung • Instrucciones de uso •
Istruzioni per l'uso • Brugsanvisning

Female • Femme •
Weiblich • Mujer •
Donna • Kvinde

Male • Homme •
Männlich • Hombre •
Uomo • Mand

empelvic®

Instructions for use • Mode d'emploi •
Gebrauchsanweisung • Instrucciones de uso •
Istruzioni per l'uso • Brugsanvisning



EMPELVIC®, THE NEW WAY TO FEEL AND TRAIN YOUR PELVIC FLOOR

Empelvic helps to give you improved awareness of your pelvic floor muscles. You will learn about the interaction between your pelvic floor, your breathing and the way your body moves. Using your cushion every day can heighten your awareness and the potential to move and control your pelvic floor.

Our instruction video on empelvic.com will help you to make the best use of your cushion. You can also visit our website where you will find more information about the pelvic floor and pelvic-floor-related conditions.

**NOTE:
INSIDE
FRONT
COVER**

SITTING COMFORTABLY ON YOUR CUSHION

- Use a chair with a flat and level seat, adjusted to your height and giving you a 90-degree bend at your knee.
- Position your Empelvic cushion as shown in the illustration.
- If you are male, sit with your anus over the anus button with the ridge behind your scrotum.
- Give yourself plenty of time to find a comfortable position. You should be able to feel the same pressure both front and back.

If it is difficult to find a good position, or if it is uncomfortable when you sit on your cushion, try the following:

- Put your knees together and spread your feet. This should give you more space under your pelvic floor.
- Place something soft under the cushion so that it adapts to fit your body, lessening the pressure on your pelvic floor.
- Remove the cushion and take a number of deep breaths, sighing as you breathe out.

Sit on the cushion again, and breathe in deeply and sigh again. Repeat this if necessary.

- Wear soft and flexible clothing.

THE EMPELVIC METHOD

The following will guide you through the various exercises, all of which can increase movement in the pelvic floor. You may not feel this movement straight away as it is very subtle and initial awareness varies from individual to individual. When you begin to feel and recognise movement in your pelvic floor during a particular exercise, keep hold of that feeling.

Carrying out these exercises can help you to become more aware of your pelvic floor and how it feels when it moves. Once you have gained this awareness you can progress to supporting and strengthening these movements, thus learning to control and steer them.

You should carry out four different exercises, each with a different purpose:

- **Exercise 1** teaches you to “let go”, allowing the pelvic floor to relax as you sigh.
- **Exercise 2** teaches you to clench actively, although still gently, as you breathe out.
- **Exercises 3 and 4** involve interaction between your pelvic floor, your stomach muscles and your breathing. You can also actively and gently strengthen the movement through clenching.

Initially, try using your Empeivic cushion for 5 min. per day.

All of the exercises can help you to achieve smoother, fluid and supple movement in your pelvic floor. Exercising every day can help you notice a difference in only a few weeks.

EXERCISE 1: SUCK AND SIGH

When you take a sharp intake of breath our pelvic floor tenses, thus giving space to allow for greater pressure in your abdomen. When you breathe out, with a sigh and a relaxed jaw, your pelvic floor relaxes.

This is the exercise you should use to massage and relax your pelvic floor. You can also use this exercise to make yourself comfortable on your cushion.

- Suck air in as if you were sucking through a straw. Do this with a degree of effort and hold on to the air. Sigh with enough force to empty your lungs gently but without breathing out with force. Avoid actively breathing out when you then exhale with a sigh. You can also make a “sighing sound” and open your mouth a little, relaxing your jaws.

Once you are used to the “suck & sigh” technique you can repeat the exercise several times and see if you can notice a response in your pelvic floor. Do you sink down a little on your cushion when you sigh? If you don't notice any difference it is not because you are doing it incorrectly, it is just because you need more time to become aware of the movement. Take your time with this exercise.

EXERCISE 2: EXHALE WITH A GENTLE SQUEEZE (CLENCH)

When you exhale, using more force, the pelvic floor will move forwards and upwards, clenching naturally. This exercise is used to tighten and close your pelvic floor. Use this exercise if you find it difficult to feel your pelvic floor, or if you want to “wake it up” and make it stronger.

- Try to support the response with an active squeeze. However, do not clench as hard as you can but use only 50% of the effort you would normally clench with. Can you feel your pelvic floor move with less difficulty when you exhale as you clench your muscles? Do you notice that a 50% squeeze feels more comfortable than a hard squeeze?

EXERCISE 3: SWIMMING STROKE

The next two exercises will teach you the interaction between your pelvic floor, your stomach muscles and your breathing.

- Place the palms of your hands together in front of your body and press them together. Move them slowly upwards, as if you are swimming. Do you notice any response in your pelvic floor as your hands and stomach muscles tense and move? Can you feel your pelvic floor relax as you lower your arms and reduce the tension in your stomach muscles?
- Try again as you breathe out while raising your arms and sigh as you lower them.

EXERCISE 4: CONSONANTS

Experiment with making different sounds as you close your eyes and focus inward.

- Start with the letter “D”. Repeat this sound loud and clear. “D.D.D.”
- Follow this with “T.T.T.T.” and, finally, “SH”, as in “hush”: “SH.SH.SH.SH.”

Perhaps you can feel a movement or “twitch” in your pelvic floor. Note the nuances in where you feel movement and just how it feels; strong, weak, movement forwards, backwards, superficial or deep.

From feeling to moving

Once you begin to notice a response in your pelvic floor to the various exercises, you know that you now have contact. You can then start to use the same exercise and focus on making your pelvic floor move more freely.

- Try to support those movements you notice in your pelvic floor with an active squeeze.

These exercises are now all about working together with your natural body movements and learning to support and strengthen them.

If you experience soreness or other discomfort, proceed with care. The cushion cannot cause any direct injury but your muscles can become a little tender when you start to use them in a different way or become newly aware of the pressure around the area. Your body will gradually become accustomed to this and the tenderness will disappear.

If you need a little more guidance, then go to our website for further information and contact details.

If any serious incident occurs with the use of this device please report to Owen Mumford Limited and your local competent authority.

Scan the QR code to see our instruction videos or find them on empelvic.com

▶▶▶▶▶ WASHING AND MATERIAL ◀◀◀◀◀

The cushion is manufactured from 100% EVA-foam. It can be cleaned using spirit or lukewarm running water and mild soap.

Dry with a cloth or leave to dry. Do not heat the cushion or you may risk burning yourself.

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Owen Mumford Ltd.,
Owen Mumford GmbH,
Freyr Life Sciences GmbH,